


IHE Integration Statement

PACS and VNA

Sectra Workstation and Sectra Core, Version 27.2, October 2025



SECTRA

Knowledge and passion

Contents

1	Introduction	1
2	About IHE	2
2.1	Integrating the Healthcare Enterprise (IHE)	2
2.2	References	2
3	IHE integration statement	3

1 Introduction

The Integrating the Healthcare Enterprise (IHE) Integration Statement was developed as a high-level statement of integration for products, which were intended to interoperate, based on IHE Integration Profiles. This document identifies the IHE Actors, IHE Integration Profiles and options Sectra AB has implemented in each product version.

The Sectra IHE Integration Statement has been created to assist the reader in determining whether and to what extent interoperability with other vendors might be supported.

The IHE Technical Framework specifies a subset of the working mechanism of the healthcare enterprise and defines their exchanges in terms of a set of coordinated transactions. The actors and transactions expressed in the IHE Technical Framework are generalizations of actual healthcare information system environments. While some of the transactions are customarily carried out by specific product groups (e.g. HIS, RIS, PACS, or modalities), the IHE Technical Framework purposely avoids relating functions or actors with such product groups. For each actor, the IHE Technical Framework identifies only those functions linked with integrating information systems. The IHE definition of an actor should not be interpreted as the absolute definition of any product that might implement that particular actor. In addition, the framework itself should not be taken as the absolute definition of healthcare information system architecture.

This IHE Integration Statement affords the reader with a sophisticated view of supported IHE integration profiles. For further exploration, additional information can be found in the DICOM and HL7 Conformance Statements of this product and in the IHE Technical Framework.

This IHE Integration Statement does not guarantee successful interoperability of this product with other vendor's products. It is the user's responsibility to thoroughly analyze the application requirements and to specify a solution that integrates the Sectra PACS or Sectra VNA with other vendor's software and systems.

This document is an addendum to the User's Guide identified on the front page. See the identified User's Guide for regulatory information.

2 About IHE

The following topics are included in this chapter:

- [Integrating the Healthcare Enterprise \(IHE\)](#)
- [References](#)

2.1 Integrating the Healthcare Enterprise (IHE)

IHE is a multi-year initiative undertaken by medical specialists, administrators, information technology professionals and manufacturers, sponsored by the Radiological Society of North America (RSNA) and the Healthcare Information and Management Systems Society (HIMSS), whose sole purpose is to improve the way computer systems in healthcare share information.

IHE is not a standards organization. Instead it promotes coordinated use of existing standards such as DICOM and HL7 to develop workflow solutions for the healthcare enterprise. Systems designed in agreement with IHE profiles communicate with one another better, are easier to implement, and facilitate the efficient access of information. Physicians, nurses, administrators and other healthcare professionals foresee a day when vital information can flow seamlessly from system to system and be readily available at the point of care. IHE is intended to make this vision a reality by improving the condition of systems integration and eliminating barriers to optimal patient care.

Though begun in the USA, the IHE is a world wide initiative. The Regional IHE Initiatives currently encompass IHE – Asia Oceania, IHE – EU (Europe) and IHE – NA (North America). The National and Regional IHE Initiatives coordinate the deployment of the IHE Technical Framework and Integration Profiles and give input to their development related to local needs and issues. They are sponsored and overseen by relevant national/regional professional associations and include membership and staff of these groups as well as local vendor representatives. These groups may rely on funding from other organizations (e.g., governmental bodies or foundations) and are responsible for securing and managing any such funds.

2.2 References

General IHE information:

- IHE home: <http://www.ihe.net>

Vendor's IHE information:

- <https://sectra.com/medical/knowledge-center/conformance-statements/>

3 IHE integration statement

Vendor	Product Name	Version	Date
Sectra AB	Sectra PACS and Sectra VNA	27.2	2025-10-10

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:

Integration Profile Implemented	Actors Implemented	Options Implemented
Cardiology Domain		
Cardiac Catheterization Workflow (CATH)	Image Manager/Image Archive	None
Echocardiography Workflow (ECHO)	Image Manager/Image Archive	None
IT Infrastructure Domain		
Patient Identifier Cross-Referencing (PIX)	Patient Identifier Cross-reference Consumer	PIX Update Notification
Consistent Time (CT)	Time Client	None
Audit Trail and Node Authentication (ATNA)	Secure Application	<ul style="list-style-type: none"> ● STX: No Secure Transport ● STX: TLS 1.2 Floor using BCP195 ● ATX: TLS Syslog ● ATX: UDP Syslog
	Audit Record Repository	<ul style="list-style-type: none"> ● ATX: TLS Syslog ● ATX: UDP Syslog
Cross Enterprise Document Sharing (XDS.b)	Document Consumer	None
	Document Source	None
	Document Administrator	Delete Metadata
Cross-Community Access (XCA)	Initiating Gateway	Affinity Domain
	Responding Gateway	None
Patient Demographics Query HL7 V3 (PDQV3)	Patient Demographics Consumer	None
Radiology Domain		

Integration Profile Implemented	Actors Implemented	Options Implemented
Scheduled Workflow (SWF)	Image Manager/Image Archive	None
	Performed Procedure Step Manager	None
Patient Information Reconciliation (PIR)	Image Manager/Image Archive	None
	Performed Procedure Step Manager	None
Consistent Presentation of Images (CPI)	Image Manager/Image Archive	None
Access to Radiology Information (ARI)	Image Manager/Image Archive	None
Key Image Note (KIN)	Evidence Creator	None
	Image Manager/Image Archive	None
Simple Image and Numeric Report (SINR)	Report Creator	None
Reporting Workflow (RWF)	Image Manager/Image Archive	None
	Performed Procedure Step Manager	None
Evidence Documents (ED)	Image Manager/Image Archive	None
Portable Data for Imaging (PDI)	Image Display	DVD Media, USB Media
	Portable Media Importer	DVD Media, USB Media
	Portable Media Creator	<ul style="list-style-type: none"> ● Web content ● DVD Media ● USB Media
NM Image (NM)	Image Manager/Image Archive	None
Teaching File and Clinical Trial Export (TCE)	Receiver	None
Cross Enterprise Document Sharing for Imaging (XDS-I.b)	Imaging Document Source	Set of DICOM Instances
	Imaging Document Consumer	None
Mammography Image (MAMMO)	Image Manager/Image Archive	None
	Image Display	None
Imaging Object Change Management (IOCM)	Change Requester	None
Cross-Community Access for Imaging (XCA-I)	Initiating Imaging Gateway	None
	Responding Imaging Gateway	None
	Imaging Document Consumer	None
Digital Breast Tomosynthesis (DBT)	Image Manager/Image Archive	<ul style="list-style-type: none"> ● For Presentation Breast Projection X-Ray Images ● For Processing Breast Projection X-Ray Images
	Image Display	For Presentation Breast Projection X-Ray Images

Integration Profile Implemented	Actors Implemented	Options Implemented
AI Results (AIR) ¹	Image Manager/Image Archive	None
Web-based Image Capture (WIC)	Image Manager	None

¹AIR supported as defined by the Rev 1.2 Trial Implementation, July 6, 2022

Compliance to Regulations and Standards

The Management System of Sectra AB [Sectra] conforms to ISO 9001, ISO 13485, ISO 27001, ISO 27017 and ISO 27018. All Sectra medical devices have obtained regulatory clearance for those markets where Sectra sells and deploys its devices, e.g. EEA, USA, Canada, Australia. For further regulatory information, please contact Sectra.

Disclaimer

Sectra shall not be liable or responsible for any use of the software in any application or manner that is not in full compliance with this manual or the Master Customer Agreement under which Sectra has made the software available to the user. This manual is intended solely as a guide and does not constitute a warranty or guarantee of functionality. All warranties and other terms and conditions governing the software are exclusively set forth in the Master Customer Agreement. Users are advised to refer to the Master Customer Agreement for comprehensive details regarding warranties and limitations of liability.

Trademarks and Patents

Sectra, Sectra Imtec and the Sectra logotype, are registered trademarks of Sectra AB. Image Exchange Portal logotype is a registered trademark of Sectra Products UK Limited. IDS7, Sectra BizTrack, Sectra CloudFlex, Sectra CTMA, Sectra DoseTrack, Sectra IEP, Sectra IMA, Sectra Image Central, Sectra Image Lab, Sectra LiteView, Sectra UniView, Sectra One, Sectra One Connect, Sectra Open Archive, Sectra VNA, Sectra UserInfluence, WISE, Sectra Upload & Store App, Sectra Upload & Store Web, Sectra Upload & Store Auto and Sectra Access Portal are trademarks of Sectra AB. Sectra Preop Online and Sectra OneScreen are trademarks of Sectra AB.

Windows is a registered trademark of Microsoft Corporation in the United States and other countries. *iPadOS*, *Mac*, *Safari*, *iPhone*, and *iPad* are registered trademarks of Apple Inc. in the United States and other countries. *Android* is a trademark of Google LLC.

All other names/products with © / ™ are registered/unregistered trademarks of the respective manufacturer.

The intellectual property of Sectra includes a number of patents, for full information refer to: <https://sectra.com/medical/about-sectra/patents/>

Copyright and Legal

© Sectra AB, Sweden, 2025

All rights are reserved. Reproduction or transmission in whole or part, in any form or by any means, electronic, mechanical or otherwise, is prohibited without written consent of the copyright owner.

Copyrights and all other proprietary rights in any software and related documentation ("Software") made available to you rest exclusively with Sectra AB. No title or ownership in the Software is conferred to you. Use of the Software is subject to the end user license conditions as are available on request.

To the maximum extent permitted by law, you shall not decompile and/or reverse engineer the software or any part thereof.

Sectra AB
Teknikringen 20
SE-58330 Linköping
Sweden
Phone: +46 13 23 52 00
E-mail: info.medical@sectra.com

SECTRA

Knowledge and passion

For other regions please visit www.sectra.com/medical/